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

## LITERATUR

1 McGuire MK, Scheyer ET, Nunn M. Evaluation of human recession defects treated with coronally advanced flaps and either enamel matrix derivative or connective tissue: comparison of clinical parameters at 10 years. J Periodontol. 2012 Nov;83(11):1353-62, Clin Oral Investig. 2016 Sep;20(7):1791-800. 2 Cairo F, Cortellini P, Pilloni A, Nieri M, Cincinelli S, Amunni F, Pagavino G, Tonetti MS. Clinical efficacy of coronally advanced flap with or without connective tissue graft for the treatment of multiple adjacent gingival recessions in the aesthetic area: a randomized controlled clinical trial. J Clin Periodontol. 2016 Oct;43(10):849-56. 3 Shirakata Y, Sculean A, Shinohara Y, Sena K, Takeuchi N, Bosshardt DD, Noguchi K. Healing of localized gingival recessions treated with a coronally advanced flap alone or combined with an enamel matrix derivative and a porcine acellular dermal matrix: a preclinical study. 4 Tonetti MS, Jepsen S; Working Group 2 of the European Workshop on Periodontology. Clinical efficacy of periodontal plastic surgery procedures: consensus report of Group 2 of the 10th European Workshop on Periodontology. J Clin Periodontol. 2014 Apr; 41 Suppl 15:S36-43. 5 Pilloni et al, Root coverage with a coronally positioned flap used in combination with enamel matrix derivative: 18 month clinical evaluation. J Periodontol. 2006; 77:2031-2039. 6 Cairo F, Nieri M, Pagliaro U. Efficacy of periodontal plastic surgery procedures in the treatment of localized facial gingival recessions. A systematic review. J Clin Periodontol. 2014 Apr;41 Suppl 15:S44-62. 7 Cairo F, Pagliaro U, Nieri M. Treatment of gingival recession with coronally advanced flap procedures: a systematic review. J Clin Periodontol. 2008 Sep;35(8 Suppl):136-62. 8 McGuire MK, Scheyer ET, Schubach P. A Prospective, Cased-Controlled Study Evaluating the use of Enamel Matrix Derivative on Human Buccal Recession Defects: A Human Histologic Examination. J Periodontol. 2016 Feb 1:1-34. 9 Jepsen S. et al. J Periodontol 2004: A Randomized Clinical Trial Comparing Enamel Matrix Derivative and Membrane Treatment of Buccal Class II Furcation Involvement in Mandibular Molars. Part I: Study Design and Results for Primary Outcomes. 10 Cieslik-Wegemund M et al, Coverage of multiple gingival recessions using the tunnel technique and collagen matrix – Mucoderm®. Poster Europerio 2015. 11 Cosgarea R, Juncor R, Lascu L, Arweiler N, Sculean A. Treatment of multiple gingival recession defects with a new collagen membrane. Poster Europerio 2015. 12 Milinkovic I, Rakaservic D, Aleksic Z, Jankovic S; Clinical application of collagen tissue matrix and enamel matrix derivative in the treatment of multiple gingival recessions. Poster Europerio 2015. 13 Horváth A, Molnár B, Gera I, Windisch P; Comparison of different approaches aimed at increasing peri-implant keratinised mucosa. Poster ITI World Symposium 2014. 14 Puisys A, Vindasiute E, Linkeviciene L, Linkevicius T. The use of acellular dermal matrix membrane for vertical soft tissue augmentation during submerged implant placement: a case series. Clin Oral Implants Res. 2015 Apr;26(4):465-70. 15 Puišys A, Žukauskas S, Kubilius R, Vindašiūtė E, Linkevicius T; Early implant placement in aesthetic area with simultaneous guided bone regeneration and soft tissue augmentation using collagen tissue matrix membrane. Clin. Oral Impl. Res. 26, Suppl. 12, 2015. 16 Schmitt CM, Matta RE, Moest T, Humann J, Gammel L, Neukam FW, Schlegel KA. Soft tissue volume alterations after connective tissue grafting at teeth: the subepithelial autologous connective tissue graft versus a porcine collagen matrix - a pre-clinical volumetric analysis. J Clin Periodontol. 2016 Jul;43(7):609-17. 17 Zafropoulos GG, Deli G, Hoffmann O, John G. Changes of the peri-implant soft tissue thickness after grafting with a collagen matrix. J Indian Soc Periodontol 2016

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© Institut Straumann AG, 2017. Alle Rechte vorbehalten.
Institut Straumann AG
### botiss biomaterials GmbH
Straumann GmbH