

# Kompatibilitätsliste für Createch

## individuelle Abutments

EXCELLENCE LINE			
Implantat	Verbindung	Plattform	Abgewinkelter Schraubenkanal (Individuelle Abutments)
ALPHA BIO®	Internal Hex	3.75	YES
AVINENT®	Internal	3.5	YES
AVINENT®	Internal	4.2	YES
AVINENT®	External	3.5	YES
AVINENT®	External	4.1	YES
AVINENT®	External	5.1	YES
BIOHORIZONS®	Internal Hex	3.0	NO
BIOHORIZONS®	Internal Hex	3.5	YES
BIOHORIZONS®	Internal Hex	4.5	YES
BIOHORIZONS®	Internal Hex	5.7	YES
BIOHORIZONS®	External Hex	4.0	YES
BIOMET 3i®	Certain™	3.4	NO
BIOMET 3i®	Certain™	4.1	YES
BIOMET 3i®	Certain™	5.0	YES
BIOMET 3i®	Certain™	6.0	YES
BIOMET 3i®	Micro mini implant	3.4	YES
BIOMET 3i®	External Hex	4.1	YES
BIOMET 3i®	External Hex	5.0	YES
BIOMET 3i®	External Hex	6.0	YES
BIOTECH®	Biotech Contact	3.3	YES
BTI®	External Hex (Tiny®)	3.5	YES
BTI®	External Hex	4.1	YES
BTI®	External Hex	5.5	YES
BTI®	Internal connection	4.1	YES
BTI®	Internal connection	5.5	YES
DENTSPLY®	Astra Tech™	Yellow	NO
DENTSPLY®	Astra Tech™	Aqua	NO
DENTSPLY®	Astra Tech™	Lilac	YES
DENTSPLY®	Xive®	Frialit™ 3.0	NO
DENTSPLY®	Xive®	Frialit™ 3.4	NO
DENTSPLY®	Xive®	Frialit™ 3.8	YES
DENTSPLY®	Xive®	Frialit™ 4.5	YES
DENTSPLY®	Xive®	Frialit™ 5.5	YES
IMPLANT DIRECT™	RePlant / Replus™	3.5	YES
IMPLANT DIRECT™	RePlant / Replus™	4.3	YES
IMPLANT DIRECT™	RePlant / Replus™	5.0	YES
IMPLANT DIRECT™	RePlant / Replus™	6.0	YES
IMPLANT DIRECT™	Legacy	3.5	YES
IMPLANT DIRECT™	Legacy	4.5	YES
KLOCKNER®	Essential®	4.5	YES
KLOCKNER®	Essential®	6.0	YES
LASAK BIONIQ®	internal Hex	4	YES
LIFECORE®	External Hex	4.1	YES
LIFECORE®	External Hex	5.0	YES
LIFECORE®	Internal Hex	3.75	YES

EXCELLENCE LINE			
Implantat	Verbindung	Plattform	Abgewinkelter Schraubenkanal (Individuelle Abutments)
LIFECORE®	Internal Hex	4.5	YES
MEGAGEN®	Exfeel External®	3.5	YES
MEGAGEN®	Exfeel External®	4.1	YES
MEGAGEN®	Exfeel External®	5.0	YES
MEGAGEN®	Ez plus™ Internal	3.5	YES
MEGAGEN®	Ez plus™ Internal	4.0	YES
MEGAGEN®	Rescue Internal®	5.0	YES
MEGAGEN®	Rescue Internal®	6.0	YES
MEGAGEN®	Anyridge® implant level	Anyridge® implant level	NO
MICRODENT®	Universal	4.10	YES
MICRODENT®	Universal	5.10	YES
MIS®	Seven	3.75	YES
MIS®	Seven	4.5	YES
MOZO GRAU®	MG InHex®	STD 3.5	YES
MOZO GRAU®	MG InHex®	MAXI 4.5	YES
NOBEL BIOCARE®	Active™	3.0	NO
NOBEL BIOCARE®	Active™	NP 3.5	NO
NOBEL BIOCARE®	Active™	RP 3.9	YES
NOBEL BIOCARE®	Brånemark System™	3.5	YES
NOBEL BIOCARE®	Brånemark System™	4.1	YES
NOBEL BIOCARE®	Brånemark System™	5.0	YES
NOBEL BIOCARE®	Replace™	3.5	NO
NOBEL BIOCARE®	Replace™	4.3	YES
NOBEL BIOCARE®	Replace™	5.0	YES
NOBEL BIOCARE®	Replace™	6.0	YES
PHIBO®	TSH	S3	YES
PHIBO®	TSH	S4	YES
PHIBO®	TSH	S5	YES
PROCLINIC®	External Hex	S3 Ø4.0	YES
PROCLINIC®	External Hex	S4 Ø4.0	YES
PROCLINIC®	External Hex	S5 Ø4.0	YES
YESGNO VINCES®	External Hex (Duo e Inttegra)	4.1	YES
S.I.N®	External Hex	4.1	YES
SOUTHERN IMPLANTS	External Hex	4.1	YES
STERNGOLD SQP®	External Hex	4.1	YES
SWEDEN & MARTINA®	External Hex (Outlink)	4.1	YES
TBR®	Hexágono externo Hex-Conic	4.0	YES
THOMMEN	PF	3.5	NO
THOMMEN	PF	4.0	NO
THOMMEN	PF	4.5	YES
THOMMEN	PF	5.0	YES
THOMMEN	PF	6.0	YES
ZIMMER®	Tapered Screw Vent & Screw Vent	3.5	YES
ZIMMER®	Tapered Screw Vent & Screw Vent	4.5	YES
ZIMMER®	Tapered Screw Vent & Screw Vent	5.7	YES

(\* All trademarks and company names are property of their respective owners.